



Installing the Killgerm Bedbug Monitor

The Bedbug Monitor is a passive trap that allows continuous monitoring for bedbug activity. It provides an attractive harbourage for bedbugs, whilst being discreet for installation in bedrooms.

The trap consists of four pieces of corrugated card, presented in such a way as to facilitate the entry of bedbugs. The white outer case is particularly appropriate for a trap of this type as 'blood spots', frequently produced by bedbugs, can easily be identified on the white surface of the outer case. It is recommended that the Bedbug Monitor is placed in areas where bedbugs are likely to be present, such as behind bed headboards, on the sides of mattresses, underneath beds, behind bedside cabinets etc.

Instructions for Use

1. The traps are pre-assembled for ease of use.
2. Using the double-sided stickers provided, place one sticker on the wider side of the trap without windows.
3. Write the installation date on the narrow side of the trap.
4. Peel the sticker to expose its sticky surface and position the trap in locations where bedbugs are likely to be present.

Killgerm Chemicals Limited, PO Box 2, Ossett, West Yorkshire, WF5 9NA.

TEL: 01924 268 400 FAX: 01924 264 757 EMAIL: sales@killgerm.com

www.killgerm.com

A member of the Killgerm Group of Companies.